

LEGIONELLA GVPC SELECTIVE SUPPLEMENT

VIAL CONTENTS (for 500 mL of medium)

Glycine	1.5 g
Vancomycin HCl	0.5 mg
Polymyxin B	40000 IU
Cycloheximide	40.0 mg

DIRECTIONS

Legionella Selective Agar

Suspend 12.5 g of Legionella BCYE Agar Base (cat. N°401582) in 500 mL of cold distilled water. Heat to boiling with agitation and sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C and add the contents of one vial Legionella BCYE α Growth Supplement reconstituted with 50 mL of sterile, warm distilled water (50°C) and the contents of one vial of Legionella GVPC Selective Supplement (cat. N° 423215) reconstituted with 10 mL of sterile distilled water. Mix well and distribute into sterile Petri dishes.

Final pH 6.9 \pm 0.1

Legionella Agar

Suspend 12.5 g of Legionella BCYE Agar Base (cat. N°401582) in 500mL of cold distilled water. Heat to boiling with agitation and sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C and add the contents of one vial Legionella BCYE α Growth Supplement (cat. N° 423210) reconstituted with 50 mL of sterile, warm distilled water (50°C). Mix well and distribute into sterile Petri dishes.

INTENDED USE

For the detection of *Legionella* spp. in environmental and clinical specimens. For the details of the procedure see the technical sheet of Legionella BCYE Agar Base (cat. N°401582)

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

Consult the material safety data sheet before the use.

PACKAGING

423215 Legionella GVPC Selective Supplement 4 vials for 500 mL of medium