

CLOSTRIDIUM DIFFICILE ANTIMICROBIC SUPPLEMENT

Vial contents (for 500 ml of medium)

D-Cycloserine.....250 mg
Cefoxitin 8 mg

DIRECTIONS

Aseptically reconstitute the contents of one vial with 5 ml of sterile distilled water. Prepare 500 ml of Clostridium Difficile Agar Base, autoclaved and cooled to 50°C. Add the supplement together with 25 ml of defibrinated sheep blood, mix well and distribute into sterile Petri dishes.

INTENDED USE

For the selective isolation of *Clostridium difficile* in clinical specimens. For the details of the procedure see the technical sheet of Clostridium Difficile Agar Base, cat. N° 401308.

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.
Consult the material safety data sheet before the use.

PACKAGING

4240006 Clostridium Difficile Antimicrobial Supplement, 10 vials, each for 500 ml of medium