

CNA ANTIMICROBIC SUPPLEMENT

VIAL CONTENTS (for 500 ml of medium)

Colistin 5.0 mg

Nalidixic Acid 7.5 mg

DIRECTIONS

Aseptically reconstitute the contents of one vial with 5 ml (or 50) of sterile distilled water. Prepare 500 ml of Columbia Agar Base or other suitable medium autoclaved and cooled to 50°C. Add the supplement together with 25 ml of defibrinated sheep blood, mix well and distribute into sterile Petri dishes.

INTENDED USE

For the selective isolation of streptococci and staphylococci. For the details of the procedure see the technical sheet of Columbia Agar Base and Columbia CNA Agar Base (cat. N° 401130-4011361).

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

Consult the material safety data sheet before the use.

PACKAGING

4240018 CNA Antimicrobial Supplement, 10 vials, each for 500 ml of medium.