

## TTC, 1% SOLUTION

**BOTTLE CONTENTS**

Triphenyltetrazolium chloride.....0.33 g

Distilled water.....30 ml

**DIRECTIONS**

The solution is ready to use.

**INTENDED USE**

For the isolation and presumptive identification of enterococci by means of the reduction of TTC to formazan, producing pink-red colonies. For the details of the procedure see the technical sheets of Azide Maltose Agar KF (cat. N° 401107)

**STORAGE**

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

**PRECAUTIONS AND WARNING**

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

**Packaging****42111801      TTC, 1% Solution,      30 ml bottle**