

ROSOLIC ACID

VIAL CONTENTS

Rosolic acid 1 g

DIRECTIONS

Transfer the contents of one vial to 100 ml volumetric flask. Add 100 ml of 0.2 N NaOH and stir to dissolve completely. Use this solution for the preparation of the Faecal Coliform Broth and faecal Coliform Agar adding 10 ml of rosolic acid solution to 1000 ml of the basal medium and boiling until complete dissolution is obtained.

INTENDED USE

Faecal Coliform Broth and Faecal Coliform Agar completed by addition of rosolic acid are used for the determination of the faecal coliforms in waters according to the membrane filtration method, with incubation at 44°C for 24 hours.

For the details of the procedure see the technical sheets of Faecal Coliform Broth/Agar, cat. N° 401486, 401487

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

Consult the material safety data sheet before the use.

PACKAGING

4211901 Rosolic Acid 1 g, 10 vials.