

## BRILLIANT GREEN 0.1% SOLUTION

**VIAL CONTENTS**

Brilliant Green 50 mg

Distilled water 50 ml

**DIRECTIONS**

Brilliant Green 0.1% Solution is ready to use

**INTENDED USE**

Brilliant Green 0.1% Solution is used to supplement Tetrathionate Broth Base (10 ml of solution/litre of base medium) and Mueller Kauffmann Tetrathionate Broth base (9.5 ml of solution/litre of base medium). For the details of the procedures see the technical sheets of Tetrathionate Broth Base (cat. N° 402125) and Mueller Kauffmann Tetrathionate Broth Base (cat. N° 401743)

**STORAGE**

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

**PRECAUTIONS AND WARNING**

The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

**PACKAGING****421505 Brilliant Green 0.1% Solution 50 ml**