

POTASSIUM TELLURITE 1% SOLUTION

Bottle contents:

Potassium Tellurite.....0.3 g

Distilled water.....30 ml

DIRECTIONS

The ready to use solutions should be used for the supplementation of the following media:
Baird Parker media, Monsur TTCA Agar Base, Giolitti Cantoni Staphylococci Broth, Vogel Johnson Agar.

INTENDED USE

For the detection of microorganisms, which reduce the potassium tellurite to tellurium, producing grey-black colonies. For the details of the procedures see the technical sheets of the above culture media.

STORAGE

Store at 2-8° - When stored as directed the supplements remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

PACKAGING**42211501****Potassium Tellurite 1% Solution,****30 ml**