

POTASSIUM TELLURITE 3.5% SOLUTION

Bottle contents:

Potassium Tellurite.....1.05 g

Distilled water.....30 ml

DIRECTIONS

The ready to use solutions should be used for the supplementation of the following media:

Baird Parker media, Monsur TTCA Agar Base, Giolitti Cantoni Staphylococci Broth, Vogel Johnson Agar.

INTENDED USE

For the detection of microorganisms, which reduce the potassium tellurite to tellurium, producing grey-black colonies. For the details of the procedures see the technical sheets of the above culture media.

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

Consult the material safety data sheet before the use.

PACKAGING

42211502 Potassium Tellurite 3.5% Solution, 30 ml