

AEROMONAS SELECTIVE SUPPLEMENT

VIAL CONTENTS (for 500 ml of medium)

Ampicillin.....5 mg

DIRECTIONS

Dissolve the contents of one vial with 5 ml of sterile distilled water, using aseptic conditions. Add to 500 ml of m-Aeromonas Selective Agar Base (Havelaar) cat. N°401019 or Aeromonas Starch DNA Agar Base, cat. N° 401020, autoclaved and cooled to 50°C. Mix well and pour into sterile Petri dishes.

INTENDED USE

For the detection of *Aeromonas* in water or other liquid samples by MF procedure, or in clinical and food samples with surface streaking technique. For the details of the procedure see the technical sheet of the two above media.

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

Consult the material safety data sheet before the use.

PACKAGING

4240012 Aeromonas Selective Supplement (Ampicillin), 10 vials, each for 500 ml of basal medium