

SALMONELLA SELECTIVE SUPPLEMENTS (FOR CHROMOGENIC MEDIUM)

VIAL CONTENTS**Salmonella Selective Supplement Vial A (vial for 500 ml)**

Emulsifying agents5,7 ml

Salmonella Selective Supplement Vial B (vial for 500 ml)

Cefsulodin..... 2,5 mg

DIRECTIONS

Prepare 500 ml of Chromogenic Salmonella Agar Base cat. N° 405350), autoclaved and cooled to 50°C; add the contents of one vial of Salmonella Selective Supplement Vial A, heat to boiling and sterilize by autoclaving at 121°C for 15 minutes. Cool to 50°C and add the contents of one vial of Salmonella Selective Supplement Vial B reconstituted with 2 ml of sterile distilled water. Mix well and pour into sterile Petri dishes.

INTENDED USE

For the isolation and identification of *Salmonella* spp For the details of the procedure see the technical sheet of Chromogenic Salmonella Agar Base (cat. N° 405350)

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

PACKAGING

4240013 Salmonella Selective Supplement 5+5 vials (each pair for 500 ml of medium)