

PP PSEUDOMONAS SUPPLEMENT

VIAL CONTENTS (for 500 ml of medium)

Penicillin G	50000UI
Pimaricin	5 mg

DIRECTIONS

Aseptically reconstitute the contents of one vial with 2 ml of sterile distilled water, mix well. Prepare 500 ml of Pseudomonas Agar Base (cat. N° 401960) autoclaved and cooled to approximately 50°C. Add the reconstituted PP Pseudomonas Supplement, mix well and pour into sterile Petri dishes.

INTENDED USE

For the selective isolation and enumeration of *Pseudomonas spp.* in milk and milk products (ISO/TS 11059). For the details of the procedure see the technical sheet of Pseudomonas Agar Base (cat. N° 401960).

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

PP Pseudomonas Supplement is classified as Xn (harmful). Consult the material safety data sheets before the use.

REFERENCE

ISO/TS 11059:2009 Milk and milk products — Method for the enumeration of *Pseudomonas spp.*

PACKAGING

4240048	PP Pseudomonas Supplement	10 vials, each for 500 ml of medium
----------------	----------------------------------	--