

PENICILLIN G ANTIMICROBIC SUPPLEMENT

Vial contents (for 100 ml of medium)

Penicillin G, sodium salt.....500 IU

DIRECTIONS

Aseptically reconstitute the contents of one vial with 2 ml of sterile distilled water. Prepare 100 ml of Sugar Free Agar Base (cat. N° 402098), autoclaved and cooled to 50°C. Add the supplement mix well and distribute into sterile Petri dishes.

INTENDED USE

For the selective enumeration of Gram negative bacilli not related to production process. For the details of the procedure see the technical sheet of Sugar Free Agar Base (cat. N° 402098).

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.
Consult the material safety data sheet before the use.

PACKAGING

4240050 Penicillin G Antimicrobial Supplement, 10 vials, each for 100 ml of medium