

## KANAMYCIN SELECTIVE SUPPLEMENT

**VIAL CONTENTS (for 500 ml of medium)**

Kanamycin Sulphate.....10 mg

**DIRECTIONS**

Dissolve the contents of one vial with 5 ml of sterile distilled water, using aseptic conditions. Add to 500 ml of Kanamycin Aesculin Azide Agar (cat. N° 401552) or to 500 ml of E.C.O.A.gar (cat. n° 401430), autoclaved and cooled to 50°C. Mix well and pour into sterile Petri dishes.

**INTENDED USE**

For the detection of enterococci in foodstuffs and other specimens. For the details of the procedure see the technical sheets of Kanamycin Aesculin Azide Agar (cat. N° 401552) or E.C.O.A.gar (401430)

**STORAGE**

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

**PRECAUTIONS AND WARNING**

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.  
Consult the material safety data sheet before the use.

**PACKAGING**

**4240055 Kanamycin Selective Supplement, 10 vials, each for 500 ml of medium**