

## POTASSIUM CHLORATE SUPPLEMENT

**Vial contents (for 250 ml of medium)**

Potassium Chlorate, 5% solution.....5 ml

**DIRECTIONS**

Prepare 250 ml of *Yersinia* ITC Broth Base (cat. N° 402270), autoclaved and cooled to 50°C. Add the Potassium Chlorate solution to the basal medium together with the contents of one vial of Ticarcillin Irganon Antimicrobial Supplement, mix well and distribute into sterile tubes.

**INTENDED USE**

For the enumeration of *Yersinia enterocolitica* in foods. For the details of the procedure see the technical sheet of *Yersinia* ITC Broth Base (cat. N° 402270).

**STORAGE**

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

**PRECAUTIONS AND WARNING**

The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

Consult the material safety data sheet before the use.

**PACKAGING**

**4240065 Potassium Chlorate Supplement, 10 vials, each for 250 ml of medium**