

## CFC PSEUDOMONAS SUPPLEMENT

### VIAL CONTENTS (for 500 ml of medium)

Cetrimide	5 mg
Fusidic Acid	5 mg
Cefalosporin	25 mg

### DIRECTIONS

Aseptically reconstitute the contents of one vial with 2 ml of sterile distilled water/ethanol (1:1), mix well. Prepare 500 ml of *Pseudomonas* Agar Base (cat. N° 401960) autoclaved and cooled to approximately 50°C. Add the reconstituted CFC Supplement, mix well and pour into sterile Petri dishes.

### INTENDED USE

For the selective isolation and enumeration of *Pseudomonas spp.* in meat and meat products (ISO 13720). For the details of the procedure see the technical sheet of *Pseudomonas* Agar Base (cat. N° 401960),

### STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

### PRECAUTIONS AND WARNING

The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

CFC *Pseudomonas* Supplement is classified as Xn (harmful). Consult the material safety data sheets before the use.

### REFERENCE

ISO 13720:2010 Meat and meat products — Enumeration of presumptive *Pseudomonas spp.*

### PACKAGING

**4240075**                      **CFC *Pseudomonas* Agar**                      **10 vials, each for 500 ml of medium**