

TWEEN® 20
TWEEN® 80**BOTTLE CONTENTS:****TWEEN 20**

Polysorbate 20 100 ml

TWEEN 80

Polysorbate 80 100 ml

DIRECTIONS

Tween 20 and Tween 80 are ready to use.

INTENDED USE

Tween 20 and Tween 80 are used as surface active agents. Tween 80 may be used for the supplementation of the following media: Lethen Broth AOAC 5 g/l, Rogosa Bios Agar 1 g/l, Rogosa Bios Broth 1 g/l.

For the details of the procedure see the technical sheets of the above media.

STORAGE

Store at 10-30° - When stored as directed the supplements remain stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

The supplements should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

PACKAGING

42120501	Tween 20	100 ml
42120502	Tween 80	100 ml