

LYSED HORSE BLOOD

BOTTLE CONTENTS

Lysed (laked) horse blood.....100 ml

INTENDED USE

Lysed horse blood should be used for the supplementation of a variety of culture media, including *Campylobacter* and *Corynebacteria* plating media.

STORAGE

Store at 2-8° - When stored as directed the supplement remains stable until the expiry date shown on the label. Do not use beyond stated expiry date.

PRECAUTIONS AND WARNING

For *in vitro* diagnostic use only. The supplement should be used only by adequately trained personnel with knowledge of microbiological techniques in the laboratory.

PACKAGING

90HLX100 **Lysed Horse Blood** **100 ml**